

Product datasheet for **GA109472**

Human SPOCK3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SPOCK3 CRISPRa kit - CRISPR gene activation of human SPARC (osteonectin), cwcv and kazal like domains proteoglycan 3
Symbol:	SPOCK3
UniProt ID:	Q9BQ16
Synonyms:	HSAJ1454; TES-3; TICN3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)