


Product datasheet for **GA109473**

Human STMN3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	STMN3 CRISPRa kit - CRISPR gene activation of human stathmin 3
Symbol:	STMN3
Components:	GA109473G1, STMN3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA109473G2, STMN3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA109473G3, STMN3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9NZ72</u>
Synonyms:	SCLIP
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)