

Product datasheet for **GA109477**

Human CDON activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	CDON CRISPRa kit - CRISPR gene activation of human cell adhesion associated, oncogene regulated
Symbol:	CDON
Components:	GA109477G1, CDON gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA109477G2, CDON gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA109477G3, CDON gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q4KMG0
Synonyms:	CDO; CDON1; HPE11; ORCAM
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)