

Product datasheet for **GA109478**

Human IMPG2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	IMPG2 CRISPRa kit - CRISPR gene activation of human interphotoreceptor matrix proteoglycan 2
Symbol:	IMPG2
Components:	GA109478G1, IMPG2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA109478G2, IMPG2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA109478G3, IMPG2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9BZV3</u>
Synonyms:	IPM200; RP56; SPACRCAN; VMD5
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)