


## Product datasheet for **GA109481**

### Human SHANK1 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	SHANK1 CRISPRa kit - CRISPR gene activation of human SH3 and multiple ankyrin repeat domains 1
Symbol:	SHANK1
Components:	GA109481G1, SHANK1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA109481G2, SHANK1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA109481G3, SHANK1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<a href="#">Q9Y566</a>
Synonyms:	SPANK-1; SSTRIP; synamon
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online »](#)