


Product datasheet for **GA109484**

Human TMED5 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TMED5 CRISPRa kit - CRISPR gene activation of human transmembrane p24 trafficking protein 5
Symbol:	TMED5
Components:	GA109484G1, TMED5 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA109484G2, TMED5 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA109484G3, TMED5 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q9Y3A6
Synonyms:	CGI-100; p24g2; p28
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)