

Product datasheet for **GA109485**

Human **SLC35B3** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	SLC35B3 CRISPRa kit - CRISPR gene activation of human solute carrier family 35 member B3
Symbol:	SLC35B3
UniProt ID:	<u>Q9H1N7</u>
Synonyms:	C6orf196; CGI-19; PAPST2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)