

Product datasheet for **GA109487**

Human TPRKB activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TPRKB CRISPRa kit - CRISPR gene activation of human TP53RK binding protein
Symbol:	TPRKB
UniProt ID:	<u>Q9Y3C4</u>
Synonyms:	CGI-121; CGI121
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)