

Product datasheet for **GA109498**

Human TMED7 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TMED7 CRISPRa kit - CRISPR gene activation of human transmembrane p24 trafficking protein 7
Symbol:	TMED7
Components:	GA109498G1, TMED7 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA109498G2, TMED7 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA109498G3, TMED7 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9Y3B3</u>
Synonyms:	CGI-109; p24g3; p24gamma3; p27
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)