

Product datasheet for **GA109502**

Human **CCDC53 (WASHC3)** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	WASHC3 CRISPRa kit - CRISPR gene activation of human WASH complex subunit 3
Symbol:	WASHC3
Components:	GA109502G1, CCDC53 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA109502G2, CCDC53 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA109502G3, CCDC53 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9Y3C0</u>
Synonyms:	CGI-116
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)