

Product datasheet for **GA109513**

Human FAM18B (TVP23B) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TVP23B CRISPRa kit - CRISPR gene activation of human trans-golgi network vesicle protein 23 homolog B
Symbol:	TVP23B
Components:	GA109513G1, FAM18B gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA109513G2, FAM18B gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA109513G3, FAM18B gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9NYZ1</u>
Synonyms:	CGI-148; FAM18B; FAM18B1; NPD008; YDR084C
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)