

Product datasheet for **GA109519**

Human SIA8C (ST8SIA3) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ST8SIA3 CRISPRa kit - CRISPR gene activation of human ST8 alpha-N-acetyl-neuraminide alpha-2,8-sialyltransferase 3
Symbol:	ST8SIA3
UniProt ID:	<u>O43173</u>
Synonyms:	SIAT8C; ST8SialIII
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)