

## Product datasheet for **GA109523**

### Human LAP3 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	LAP3 CRISPRa kit - CRISPR gene activation of human leucine aminopeptidase 3
Symbol:	LAP3
UniProt ID:	<a href="#">P28838</a>
Synonyms:	HEL-S-106; LAP; LAPEP; PEPS
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)