

## Product datasheet for GA109551

### Human TNNI3K activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	TNNI3K CRISPRa kit - CRISPR gene activation of human TNNI3 interacting kinase
Symbol:	TNNI3K
UniProt ID:	<a href="#">Q59H18</a>
Synonyms:	CARK; CCDD
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)