

Product datasheet for GA109555

Human SEPSECS activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SEPSECS CRISPRa kit - CRISPR gene activation of human Sep (O-phosphoserine) tRNA:Sec (selenocysteine) tRNA synthase
Symbol:	SEPSECS
UniProt ID:	Q9HD40
Synonyms:	LP; MGC161491; Sep (O-phosphoserine) tRNA:Sec (selenocysteine) tRNA synthase; SLA; SLA/LP; soluble liver antigen/liver pancreas antigen
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)