

Product datasheet for **GA109556**

Human SIDT2 activation kit by CRISPRa

Product data:

| | |
|---------------|--|
| Product Type: | CRISPRa Kits |
| Description: | SIDT2 CRISPRa kit - CRISPR gene activation of human SID1 transmembrane family member 2 |
| Symbol: | SIDT2 |
| UniProt ID: | Q8NBJ9 |
| Synonyms: | CGI-40 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)