

Product datasheet for **GA109558**

Human Adiponectin Receptor 1 (ADIPOR1) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ADIPOR1 CRISPRa kit - CRISPR gene activation of human adiponectin receptor 1
Symbol:	ADIPOR1
UniProt ID:	Q96A54
Synonyms:	ACDCR1; CGI-45; CGI45; FLJ25385; FLJ42464; PAQR1; TESBP1A
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)