

Product datasheet for **GA109586**

Human IER3IP1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	IER3IP1 CRISPRa kit - CRISPR gene activation of human immediate early response 3 interacting protein 1
Symbol:	IER3IP1
UniProt ID:	Q9Y5U9
Synonyms:	HSPC039; MEDS; PRO2309
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)