

Product datasheet for **GA109595**

Human KCTD3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	KCTD3 CRISPRa kit - CRISPR gene activation of human potassium channel tetramerization domain containing 3
Symbol:	KCTD3
Components:	GA109595G1, KCTD3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA109595G2, KCTD3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA109595G3, KCTD3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9Y597</u>
Synonyms:	NY-REN-45
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)