

Product datasheet for **GA109598**

Human PTD016 (RNFT1) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	RNFT1 CRISPRa kit - CRISPR gene activation of human ring finger protein, transmembrane 1
Symbol:	RNFT1
Components:	GA109598G1, PTD016 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA109598G2, PTD016 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA109598G3, PTD016 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q5M7Z0</u>
Synonyms:	PTD016
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)