

Product datasheet for **GA109607**

Human C5orf45 (MRNIP) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	MRNIP CRISPRa kit - CRISPR gene activation of human MRN complex interacting protein
Symbol:	MRNIP
UniProt ID:	<u>Q6NTE8</u>
Synonyms:	DKFZp686L2452; MGC65027; MGC78537
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)