

Product datasheet for **GA109610**

Human C1orf33 (MRT04) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	MRT04 CRISPRa kit - CRISPR gene activation of human MRT4 homolog, ribosome maturation factor
Symbol:	MRT04
Components:	GA109610G1, C1orf33 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA109610G2, C1orf33 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA109610G3, C1orf33 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9UKD2</u>
Synonyms:	C1orf33; dj657E11.4; MRT4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)