

Product datasheet for **GA109632**

Human GPN3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	GPN3 CRISPRa kit – CRISPR gene activation of human GPN-loop GTPase 3
Symbol:	GPN3
Components:	GA109632G1, GPN3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA109632G2, GPN3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA109632G3, GPN3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9UHW5</u>
Synonyms:	ATPBD1C
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

