

Product datasheet for **GA109636**

Human SS18L2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SS18L2 CRISPRa kit - CRISPR gene activation of human SS18 like 2
Symbol:	SS18L2
Components:	GA109636G1, SS18L2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA109636G2, SS18L2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA109636G3, SS18L2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9UHA2</u>
Synonyms:	KIAA-iso
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)