

Product datasheet for **GA109639**

Human ZNF 639 (ZNF639) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ZNF639 CRISPRa kit - CRISPR gene activation of human zinc finger protein 639
Symbol:	ZNF639
Components:	GA109639G1, ZNF 639 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA109639G2, ZNF 639 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA109639G3, ZNF 639 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9UID6</u>
Synonyms:	ANC-2H01; ANC_2H01; ZASC1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)