

Product datasheet for **GA109660**

Human PIGP activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PIGP CRISPRa kit - CRISPR gene activation of human phosphatidylinositol glycan anchor biosynthesis class P
Symbol:	PIGP
UniProt ID:	P57054
Synonyms:	DCRC; DCRC-S; DSRC; PIG-P; PIGP
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)