

Product datasheet for **GA109666**

Human TMEM85 (EMC4) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	EMC4 CRISPRa kit – CRISPR gene activation of human ER membrane protein complex subunit 4
Symbol:	EMC4
Components:	GA109666G1, TMEM85 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA109666G2, TMEM85 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA109666G3, TMEM85 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q5J8M3</u>
Synonyms:	PIG17; TMEM85
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

