

Product datasheet for **GA109669**

Human ANKRD39 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ANKRD39 CRISPRa kit - CRISPR gene activation of human ankyrin repeat domain 39
Symbol:	ANKRD39
Components:	GA109669G1, ANKRD39 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA109669G2, ANKRD39 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA109669G3, ANKRD39 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q53RE8</u>
Synonyms:	MGC41816
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)