

## Product datasheet for **GA109677**

### Human NT5C3 (NT5C3A) activation kit by CRISPRa

#### Product data:

|               |  |
|---------------|--|
| Product Type: | CRISPRa Kits   |
| Description:  | NT5C3A CRISPRa kit - CRISPR gene activation of human 5'-nucleotidase, cytosolic IIIA |
| Symbol:       | NT5C3A   |
| UniProt ID:   | <u><a href="#">Q9H0P0</a></u>  |
| Synonyms:     | cN-III; hUMP1; NT5C3; P5 N-1; P5N-1; p36; PN-I; POMP; PSN1; UMPH; UMPH1              |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)              |

#### Product images:



[View online >](#)