

Product datasheet for **GA109692**

Human TFDP3 activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | TFDP3 CRISPRa kit - CRISPR gene activation of human transcription factor Dp family member 3 |
| Symbol: | TFDP3 |
| UniProt ID: | Q5H9I0 |
| Synonyms: | CT30; DP4; HCA661 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)