

Product datasheet for **GA109693**

Human UBAP1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	UBAP1 CRISPRa kit - CRISPR gene activation of human ubiquitin associated protein 1
Symbol:	UBAP1
Components:	GA109693G1, UBAP1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA109693G2, UBAP1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA109693G3, UBAP1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9NZ09</u>
Synonyms:	NAG20; UAP; UBAP; UBAP-1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)