

Product datasheet for **GA109717**

Human THEG activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | THEG CRISPRa kit - CRISPR gene activation of human theg spermatid protein |
| Symbol: | THEG |
| UniProt ID: | Q9P2T0 |
| Synonyms: | CT56; THEG1 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)