

Product datasheet for **GA109724**

Human KCNK9 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	KCNK9 CRISPRa kit - CRISPR gene activation of human potassium two pore domain channel subfamily K member 9
Symbol:	KCNK9
UniProt ID:	Q9NPC2
Synonyms:	K2p9.1; KT3.2; TASK-3; TASK3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)