


Product datasheet for **GA109731**

Human **SLC25A37** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	SLC25A37 CRISPRa kit - CRISPR gene activation of human solute carrier family 25 member 37
Symbol:	SLC25A37
Components:	GA109731G1, SLC25A37 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA109731G2, SLC25A37 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA109731G3, SLC25A37 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9NYZ2</u>
Synonyms:	HT015; MFRN; MFRN1; MSC; MSCP; PRO1278; PRO1584; PRO2217
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)