

Product datasheet for **GA109733**

Human TXNDC3 (NME8) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	NME8 CRISPRa kit - CRISPR gene activation of human NME/NM23 family member 8
Symbol:	NME8
Components:	GA109733G1, TXNDC3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA109733G2, TXNDC3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA109733G3, TXNDC3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8N427</u>
Synonyms:	CILD6; HEL-S-99; NM23-H8; sptrx-2; SPTRX2; TXNDC3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)