

Product datasheet for **GA109737**

Human MRPL35 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	MRPL35 CRISPRa kit - CRISPR gene activation of human mitochondrial ribosomal protein L35
Symbol:	MRPL35
Components:	GA109737G1, MRPL35 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA109737G2, MRPL35 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA109737G3, MRPL35 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9NZE8</u>
Synonyms:	L35mt; MRP-L35
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)