


Product datasheet for **GA109738**

Human **RSRC1** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	RSRC1 CRISPRa kit - CRISPR gene activation of human arginine and serine rich coiled-coil 1
Symbol:	RSRC1
Components:	GA109738G1, RSRC1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA109738G2, RSRC1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA109738G3, RSRC1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q96IZ7</u>
Synonyms:	BM-011; SFRS21; SRrp53
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)