

Product datasheet for **GA109739**

Human MEX3C activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	MEX3C CRISPRa kit - CRISPR gene activation of human mex-3 RNA binding family member C
Symbol:	MEX3C
UniProt ID:	<u>Q5U5Q3</u>
Synonyms:	BM-013; MEX-3C; RKHD2; RNF194
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)