

## Product datasheet for **GA109741**

### Human WAC activation kit by CRISPRa

#### Product data:

|               |   |
|---------------|---|
| Product Type: | CRISPRa Kits  |
| Description:  | WAC CRISPRa kit - CRISPR gene activation of human WW domain containing adaptor with coiled-coil |
| Symbol:       | WAC   |
| UniProt ID:   | <u><a href="#">Q9BTA9</a></u>   |
| Synonyms:     | BM-016; DESSH; PRO1741; Wwp4  |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                         |

#### Product images:



[View online >](#)