

## Product datasheet for **GA109746**

### Human TWEAKR (TNFRSF12A) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	TNFRSF12A CRISPRa kit - CRISPR gene activation of human TNF receptor superfamily member 12A
Symbol:	TNFRSF12A
UniProt ID:	<a href="#">Q9NP84</a>
Synonyms:	CD266; FN14; TWEAKR
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)