

Product datasheet for **GA109749**

Human PRR16 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PRR16 CRISPRa kit - CRISPR gene activation of human proline rich 16
Symbol:	PRR16
Components:	GA109749G1, PRR16 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA109749G2, PRR16 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA109749G3, PRR16 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q569H4
Synonyms:	DSC54; LARGEN
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)