

## Product datasheet for **GA109757**

### Human **TAOK3** activation kit by **CRISPRa**

#### Product data:

Product Type:	CRISPRa Kits
Description:	TAOK3 CRISPRa kit - CRISPR gene activation of human TAO kinase 3
Symbol:	TAOK3
UniProt ID:	<a href="#">Q9H2K8</a>
Synonyms:	DPK; JIK; MAP3K18
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)