

## Product datasheet for **GA109770**

### Human HSPC014 (POMP) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	POMP CRISPRa kit - CRISPR gene activation of human proteasome maturation protein
Symbol:	POMP
UniProt ID:	<a href="#">Q9Y244</a>
Synonyms:	C13orf12; HSPC014; PNAS-110; UMP1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)