

Product datasheet for **GA109782**

Human IF3EI (EIF3L) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	EIF3L CRISPRa kit - CRISPR gene activation of human eukaryotic translation initiation factor 3 subunit L
Symbol:	EIF3L
Components:	GA109782G1, IF3EI gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA109782G2, IF3EI gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA109782G3, IF3EI gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9Y262</u>
Synonyms:	EIF3EIP; EIF3S6IP; EIF3S11; HSPC021; HSPC025; MSTP005
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)