

Product datasheet for **GA109783**

Human NIP7 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	NIP7 CRISPRa kit - CRISPR gene activation of human nucleolar pre-rRNA processing protein NIP7
Symbol:	NIP7
Components:	GA109783G1, NIP7 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA109783G2, NIP7 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA109783G3, NIP7 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9Y221</u>
Synonyms:	CGI-37; HSPC031; KD93
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)