


Product datasheet for **GA109798**

Human ZNF107 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ZNF107 CRISPRa kit - CRISPR gene activation of human zinc finger protein 107
Symbol:	ZNF107
Components:	GA109798G1, ZNF107 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA109798G2, ZNF107 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA109798G3, ZNF107 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9UII5</u>
Synonyms:	smap-7; Y8; ZFD25; ZNF588
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)