

Product datasheet for **GA109852**

Human C20orf111 (OSER1) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	OSER1 CRISPRa kit - CRISPR gene activation of human oxidative stress responsive serine rich 1
Symbol:	OSER1
Components:	GA109852G1, C20orf111 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA109852G2, C20orf111 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA109852G3, C20orf111 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q9NX31
Synonyms:	C20orf111; dj1183l21.1; HSPC207; Osr1; Perit1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)