

Product datasheet for **GA109873**

Human NT5DC3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	NT5DC3 CRISPRa kit – CRISPR gene activation of human 5'-nucleotidase domain containing 3
Symbol:	NT5DC3
Components:	GA109873G1, NT5DC3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GEI00074) GA109873G2, NT5DC3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GEI00074) GA109873G3, NT5DC3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GEI00074) 1 CRISPRa-Enhancer vector, SKU GEI00056 1 CRISPRa scramble vector, SKU GEI00077
UniProt ID:	<u>Q86UY8</u>
Synonyms:	TUI2B1-TY; TUI2BITY
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

