

Product datasheet for **GA109875**

Human IL23 (IL23A) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	IL23A CRISPRa kit - CRISPR gene activation of human interleukin 23 subunit alpha
Symbol:	IL23A
UniProt ID:	<u>Q9NPF7</u>
Synonyms:	IL-23; IL-23A; IL23P19; P19; SGRF
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)