

Product datasheet for **GA109878**

Human ARM CX3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ARM CX3 CRISPRa kit - CRISPR gene activation of human armadillo repeat containing X-linked 3
Symbol:	ARM CX3
Components:	GA109878G1, ARM CX3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA109878G2, ARM CX3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA109878G3, ARM CX3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9UH62</u>
Synonyms:	ALEX3; dj545K15.2; GASP6
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)